

**From the original research**

# **TECH BIOSEALER**

**the endodontic cement inducing  
formation of apatite**



**bioactivity and biocompatibility**

**isasan**

# TECH BIOSEALER

✓ Biocompatibility

✓ Bioactivity with apatite formation

✓ Antibacterial activity



✓ Harden in damp conditions and in presence of organic fluids


✓ Marginal adjustment

✓ Dimensional stability (phyllosilicate patented)


✓ Wear resistance

✓ Adequate expansion


ENDO




CAPPING



APEX



ROOT END




## BIOCOMPATIBILITY

The SEM observations show that cells are able to adhere the surface and proliferate colonizing the material (osteoconductive activity of the cements).

*Gandolfi MG et al. (J Endod 2008, JBMR 2008, J Endod 2010)*

## ANTIBACTERIAL activity



The cement increases the pH of the fluids till 12, playing then antibacterial activity. The pH remains high for several days.


*Gandolfi MG et al. (Clin Oral Invest 2009)*



Inhibition area of bacterial growth (*Enterococcus faecalis*, *Streptococcus mutans*, *Pseudomonas aeruginosa*) of Tech BIOSEALER ENDO cement.


*(Sambri V 2009 - per gentile concessione Ist. Microbiologia Università Bologna)*

# BIOACTIVITY



The surface analysis by micro-Raman spectroscopy shows the formation of an apatite film already after 7 days.

Taddei P (da Gandolfi MG et al, Clin Oral Invest 2009)



**Harden in wet conditions and with presence of organic fluids**





BIOSEALER will harden also in presence of blood, plasma, saliva and dentinal fluid.

(courtesy of Chersoni S)



# OSTEOINDUCTION

Contact with organic fluids (blood, plasma, dentinal fluid, saliva) starts the formation of apatite on the surface of Tech BIOSEALER cements. The analysis by SEM-EDX show a layer of calcium phosphate (apatite) on the surface of cements and a calcium silicated structure inside.

Gandolfi MG et al. (Clin Oral Invest 2009, JABB 2009)



## BIOACTIVITY




### Basic research.

Tech BIOSEALER cements produce spherulites of apatite (bioactivity) on the surface already after 24 hours contact with organic fluids. The apatite is the mineral that forms and builds dentin and bone.


*Gandolfi MG et al. (J Endod 2010)*

## Dimensional stability and marginal adaptation


### Tech Biosealer cement with phyllosilicate



### Tech Biosealer cement with phyllosilicate



### Cement without phyllosilicate



Tech BIOSEALER cements have dimensional stability and marginal adaptation because of their innovative formulation and the presence of a phyllosilicate patented. The electron microscopy observations show the absence of porosity and marginal gaps at the interface dentin/cement.

*Gandolfi MG et al. (J Endod 2007).*



### Cement for endodontic treatments.

- ✓ Perfect biocompatibility
- ✓ Hardens in damp or wet channels (with open apex)
- ✓ Excellent combination with warm gutta-percha
- ✓ Volumetric expansion of 0.35 to 12.50%
- ✓ Seal of the porosity thanks to the formation of apatite
- ✓ High antibacterial activity resulting from the release of calcium hydroxide
- ✓ Excellent result of the apical seal after 12 months
- ✓ Excellent RX-opacity

#### Packaging:

blister pack of 10 capsules of 270 mg each.  
1 bottle 5 ml of liquid DPBS



### Cement for retrograde fillings

- ✓ Perfect biocompatibility and bioactivity
- ✓ Immediate marginal adaptation
- ✓ Reduced hardening time after application in rear-cavity
- ✓ Hardens in any condition of humidity and also in presence of blood
- ✓ Perfect marginal stable seal
- ✓ Formation of a bioactive apatite surface
- ✓ High antibacterial activity and RX-opacity

#### Packaging:

blister pack of 5 capsules of 270 mg each.  
1 bottle 5 ml of liquid DPBS



### Cement for apexification

- ✓ Excellent biocompatibility
- ✓ High release of calcium hydroxide
- ✓ Hardens in damp conditions due to open apex
- ✓ Hardens also mixed with blood or other body fluids
- ✓ Induces the formation of bone tissue (apatite) at apical level
- ✓ It obtains a definite apical seal

#### Packaging:

blister pack of 5 capsules of 270 mg each.  
1 bottle 5 ml of liquid DPBS



### Cement for direct capping

- ✓ Hardens quickly also in presence of blood and humidity
- ✓ Immediate formation and release of calcium hydroxide
- ✓ Formation of secondary dentine in 3-4 weeks
- ✓ Perfect biocompatibility with exposed pulp tissue
- ✓ Stimulates pulp stem cells

#### Packaging:

blister pack of 5 capsules of 270 mg each.  
1 bottle 5 ml of liquid DPBS

- [1] H.K. Chang, I. Islam, E.T. Koh, J. Endod. 31 (2006) 665.
- [2] W.P. Saunders, J. Endod. 34 (2008) 660.
- [3] D.T. Holden, S.A. Schwartz, T.C. Kirkpatrick, W.G. Schindler, J. Endod. 34 (2008) 812.
- [4] R. Pace, V. Giuliani, G. Pagavino, J. Endod. 34 (2008) 1130.
- [5] M.G. Gandolfi, S. Sauro, F. Mannocci, S. Zanna, M. Capoferri, C. Prati, R. Mongiorgi, J. Endod. 33 (2007) 1082.
- [6] M.G. Gandolfi, F. Perut, G. Ciapetti, R. Mongiorgi, C. Prati, J. Endod. 34 (2008) 39.
- [7] M.G. Gandolfi, S. Farascioni, D.H. Pashley, G. Gasparotto, C. Prati, J. Dent. 36 (2008) 565.
- [8] P.J.C. Mitchell, T.R. Pitt Ford, M. Torabinejad, F. McDonald, Biomaterials 20 (1999) 167.
- [9] E.T. Koh, F. McDonald, T.R. Pitt Ford, M. Torabinejad, J. Endod. 24 (1998) 24.
- [10] A.M. Belu, D.J. Graham, D.G. Castner, Biomaterials 24 (2003) 3635.
- [11] H.B. Lu, C.T. Campbell, D.J. Graham, B.D. Ratner, Anal. Chem. 72 (2000) 2886.
- [12] P. Malmberg, U. Bexell, C. Eriksson 1, H. Nygren, K. Richter, Rapid Commun. Mass. Spectrom. 21 (5) (2007) 745.
- [13] C. Eriksson, P. Malmberg, H. Nygren, Rapid Commun. Mass Spectrom. 22 (7) (2008) 943.
- [14] D. Touboul, F. Kollmer, E. Niehuis, A. Brunelle, O. Lapr evote, J. Am. Soc. Mass Spectrom. 16 (10) (2005) 1608.
- [15] M.S. Wagner, B.J. Tyler, D.G. Castner, Anal. Chem. 74 (2002) 1824.
- [16] D.J. Graham, M.S. Wagner, D.G. Castner, Appl. Surf. Sci. 252 (2006) 6860.
- [17] C. Eriksson, K. B rner, H. Nygren, K. Ohlson, U. Bexell, N. Billerdahl, M. Johansson, Appl. Surf. Sci. 252 (2006) 6757.
- [18] N.K. Sarkar, R. Caceido, P. Ritwik, R. Moiseyeva, I. Kawashima, J. Endod. 31 (2005) 97.
- [19] P. Taddei, A. Tinti, M.G. Gandolfi, P. Rossi, C. Prati, J. Mol. Struct. 924–926 (2009).
- [20] Mohn 2010, Gandolfi OOOE 2009, Gandolfi Clin Oral Invest 2009, Gandolfi JAAB (J Applied Biomaterials Biomechanics) 2009, Gandolfi J Endodontics 2010
- [21] Taddei P, Tinti A, Gandolfi MG, Rossi PL, Prati C. Ageing of calcium silicate cements for endodontic use in simulated body fluids: a micro-Raman study. J Raman Spectrosc 2009; 40: 1858-1866.
- [22] Gandolfi MG, Ciapetti G, Perut F, Taddei P, Modena E, Rossi PL, Prati C. Biomimetic calcium-silicate cements aged in simulated body solutions. Osteoblasts response and analyses of apatite coating. J Applied Biomaterials and Biomechanics 2009; 7: 160-170.
- [23] Gandolfi Maria Giovanna, Iacono Francesco, Agee Kelly, Siboni Francesco, Tay Franklin, Pashley David Herbert, Prati Carlo. Setting time and expansion in different soaking media of experimental accelerated calcium-silicate cements and ProRoot MTA. OOOE 2009; 108: e39-e45.
- [24] Gandolfi Maria Giovanna, Taddei Paola, Tinti Anna, Dorigo De Stefano Elettra, Rossi Piermaria Luigi, Prati Carlo. Kinetics of apatite formation on a calcium-silicate cement for root-end filling during ageing in physiological-like phosphate solutions. Clin Oral Invest 2009, in press
- [25] Gandolfi MG, Van Landuyt K, Taddei P, Modena E, Van Meerbeek B, Prati C. ESEM-EDX and Raman techniques to study MTA calcium-silicate cements in wet conditions and in real-time. Joe 2010; 36: 851-857.
- [26] Torrisi A, Torrisi V, Tuccitto N, Gandolfi MG, Prati C, Licciardello A. ToF SIMS images and spectra of biomimetic calcium-silicate based cements after storage in solutions simulated the human biological fluid effects. Int J Mass Spectrosc 2010; 289: 150-161.
- [27] Mohn D, Zehnder M, Imfeld T, Srark WJ. International Endodontic Journal 2010; 43: 210-217.
- [28] Gandolfi MG, Fiorentini E, Siboni F, Devescovi V, Colin A, Marchionni S, Ciapetti G, Rossi PL, Prati C. Utilizzo di cellule staminali pulpari su scaffold a base di cementi calcio-silicati. Odontoiatria Rivista degli Amici di Brugg 2009; 1:9-12



**TECH BIOSEALER** is a trademark ISASAN

**isasan**

Via Bellini, 17 - 22070 Rovello Porro (CO) tel. 02 96 75 41 79 - fax 02 96 75 41 90 [www.isasan.com](http://www.isasan.com) - [info@isasan.com](mailto:info@isasan.com)

distributore esclusivo

**TECH**

**(FGM)**

saniswiss

**MRC**  
MYOFUNCTIONAL  
RESEARCH CO.

**StickTech**

**paki**<sub>43</sub>